



**No Disclosures**


Introduction


**POCT in Massive  
Bleedings**


Conventional  
Assays

New Assays


## Introduction









## Conventional Assays




## Pitfall 2 Analyser





## Conventional Assays

### Pitfall 2 Analyzer




#### Errors and warnings

No errors

Measured	s	Replicate 1	Replicate 2	Mean
	25.8		25.8	
	25.8		25.8	

## Conventional Assays


### Pitfall 3 Fibrinogen after haemorrhage



Fries et al., Br J Anaesth 2010

## Conventional Assays

### Pitfall 3 Fibrinogen after haemorrhage



Fries et al., Br J Anaesth 2010

## Conventional Assays


**Fibrinogen levels during trauma hemorrhage, response to replacement therapy, and association with patients outcomes**  
 Rourhe C, Curry N, Khan S, Taylor R, Raza I, Davenport R, Stanworth S, Brohi K  
*J Throm Haemost 2012*

517/555 enrolled patients


Exclusion criteria: < 2000 mL of fluid before hospital arrival

ACT defined on Rotem results

### Pitfall 3 Fibrinogen during trauma haemorrhage



## Conventional Assays



### Pitfall 3 Analyzer

#### Errors and warnings

Replicate 1		
Measured (s)		
CE	5071	(Coag) First derivative peak not found
		Unit 1 (s)
RE	5100	Measured result failed

#### Validation rules


Rerun: modified fibrinogen, FibLow...

**Coagulation results obtained after at least 45 min**

## New Assays




### Thrombin Generation Assay (TGA)



## Conventional Assays

### Pitfall 2 Steps analysed with CCT



## Conventional Assays

### TGA and Massive Bleeding

Little data in human  
No POCT  
Requires trained staff  
Independent of fibrinogen level  
No transfusion chart described  
In vitro studies  
Mainly studies on animal trauma models

## Thrombin Generation Assay (TGA)

### TGA and Resuscitation Fluids

Altered TGA with HES but not with other fluids  
Correction with PCC or cryoprecipitate  
Weak improvement with Novoseven  
No correction with fibrinogen

### TGA and Acidosis



Pig trauma model  
Induced acidosis  
Reduction of the ETP and thrombin peak


*Caballo et coll, Blodd Transfu 2012  
Darlington et coll, J Trauma 2011*

Rotem

New Assays

## Thromboelastography/ Thromboelastometry





## New Assays

### **Admission rapid thrombelastography can replace conventional coagulation tests in the emergency department: experience with 1974 consecutive trauma patients.**

Holcomb J, Minei K, Scerbo M, Radwan Z, Wade C, Kozar R, Gill B, Albarado R, McNutt M, Khan S, Adams P, McCarthy J, Cotton B  
*Ann Surg* 2012

#### **Patients**

Memorial Hermann Hospital, Houston, Trauma Center

Deep trauma, exclusion of third-degree burns

>18 y

September 2009 – February 2011

## Rotem

#### **Analysis**

Immediately at the admission

PT, aPTT, Fibrinogen      }  
Platelets                    } CCT (Conventional Coagulation Tests)  
r-TEG

#### **Outcomes**

Correlations r-TEG-CCT (<0.3, 0.3-0.7, >0.7)

Correlations r-TEG-transfusions

Correlations CCT-transfusions

## New Assays

### Admission rapid thrombelastography can replace conventional coagulation tests in the emergency department: experience with 1974 consecutive trauma patients.

Holcomb J, Minei K, Scerbo M, Radwan Z, Wade C, Kozar R, Gill B, Albarado R, McNutt M, Khan S, Adams P, McCarthy J, Cotton B  
*Ann Surg 2012*

#### Résultats

1908 trauma patients included

#### Correlation of r-TEG values with CCTs (n=1974)

## Rotem

	<u>PT</u>	<u>aPTT</u>	<u>INR</u>	<u>Platelet Count</u>	<u>Fibrinogen</u>
ACT, sec	0.35	0.47	0.52	-0.15	-0.17
r-value, min	0.24	0.32	0.37	-0.14	-0.17
k-time, min	0.21	0.44	0.34	-0.25	-0.32
α-angle, degree	-0.23	-0.41	-0.33	0.34	0.53
MA, mm	-0.22	-0.35	-0.27	0.42	0.63
G-value	-0.02	-0.03	-0.03	-0.01	0.01

*ACT, r-value:* deficit in factors, hemodilution

*k-time:* deficit in fibrinogen, (deficit in factors)

*α-angle, degree:* deficit in fibrinogen, deficit in platelets

*MA:* deficit in fibrinogen, deficit in platelets

*G-value:* clots' strength

## New Assays

### **Admission rapid thrombelastography can replace conventional coagulation tests in the emergency department: experience with 1974 consecutive trauma patients.**

Holcomb J, Minei K, Scerbo M, Radwan Z, Wade C, Kozar R, Gill B, Albarado R, McNutt M, Khan S, Adams P, McCarthy J, Cotton B  
*Ann Surg 2012*

*Multiple regression analysis*

Prediction of pRBC transfusion in the first 6 hours of admission:

Predicted: r-TEG, PT, aPTT, INR

Non predicted: fibrinogen, platelets

Prediction of FFP and platelets transfusion

Predicted: r-TEG, PT, aPTT, INR, platelets number

Non predicted: fibrinogen

Massive transfusion

Predicted: r-TEG and CCT ( $\alpha$  angle is the more predicted)

Mortality at 6h and 30d

Predicted: r-TEG, aPTT

Non predicted: PT, INR, platelets number, fibrinogen

### **Conclusions**

r-TEG is better than the 5 CCT in the management of MB

Early detection of patient with acute needs of transfusion

r-TEG at the admission could (should) replace CCT

*European Recommendations: grade 2c (2010) versus grade 1c (2013)*

## New Assays

### Rotem

#### **Advantages**

Extensive coagulation exploration  
Coagulation cascade  
Fibrino-formation  
Fibrinolysis

Whole blood  
Short TAT  
No centrifugation  
Hematocrit  
Platelets

POCT





Real-time analysis, no buffering  
Transfusion protocols  
Recommended in the European guidelines

#### **Disadvantages**

Not sensitive to platelets function  
Inhibitors of platelet function  
Must be used with platelet function analyser (Verify-Now...)  
Moderately sensitive to anticoagulants (heparin, ...)  
Performed at 37° C  
Recalcification  
High concentration of TF  
No information about the sub-endothelium functionnality

## New Assays

### Multiplate



**Solomon et coll, Thromb Haemos 2011**

Retrospective study on 163 traumas (ISS>18, 12,3% death)


ADP, COL, TRAP

ADP and TRAP: abnormal values in the non survival group

Platelet dysfunction is a sign of ACT

## New Assays

### Thromboelastography/ Thromboelastometry



## POCT in Massive Bleedings

### Conventional Assays

Specific assays  
Slow  
No transfusion protocols  
Not predictive enough

Accurate  
Validated  
Cheap

pRBC: 116 euros  
FFP: 90 euros  
Platelets: 413 euros  
Fibrinogen: 350 euros

### New Assays

Global assays  
Fast  
Real-time results  
Transfusion protocols  
Predictive of blood requirement

Validation  
Variability  
More expensive